

2020

Guide for CIA

Proceedings

David Pérez González

22.12.2020

Contents

1. Aim of Proceeding	3
2. Introduction	3
3. Obtaining Data	3
4. Loading Data	4
5. Analysing Data	5
5.1. Peak Finder	5
5.2. First Part of CIA	5
5.2.1. Assignment of Label to Crystals	6
5.2.2. First Identification	7
5.2.3. Classification of Crystal	7
5.3. Second Part of CIA	11
6. Cluster Identification	13
A. Settings Measurement.ini	15

1. Aim of Proceeding

This document contains the guide for explaining how to use the Crystal Identification Algorithm (CIA) in Python. The document itself conforms to its own specifications. This guide is done to avoid possible issues with the use of this feature and to take care of the constraints that one needs to take into account.

2. Introduction

The following guide is based on the beta version of the CIA and on the version 2.7 of Python.

The goal of this algorithm is to be able to calibrate the PET-detectors by identifying the different detected crystals and with that, obtain DOI information which shall improve the energy and time resolution of our detector.

Four flood maps are generated from four different center-of-gravity (COG) algorithms out of the raw data. Each of the flood maps is split into smaller regions of interest (ROIs) where the relative coordinates of the clusters are concentrated. Then, a peak-finder routine is run through each ROI to look for the clusters belonging to each crystal and the center-peak of the ROI is labelled using a LUT.

With this label, the rest of the peaks from the ROI are as well identified going towards the edges. A double-check is applied later with all the ROIs together from the same flood map. Medians for the rows and columns are calculated and based on them, the peaks are confirmed as valid.

With all the labels from the four COGs we expect to have the complete spectrum of crystals identified. A LUT is created based on this information. All the clusters are then assigned to the crystals using the LUT. In principle, a homogeneous distribution of the clusters is expected, but some minor values from one of the COG algorithms are found both in this distribution and its corresponding flood map.

3. Obtaining Data

To obtain the positions of the clusters from the four different COGs, the hitAnalysis program (from mecfor satumanalysis) has to be run again with the corresponding "measurement.txt" file and in the folder where the data was recorded. All the measurements are under the following folder:

- nas/pet-scratch/Measurements/Hypmed/

The command itself is:

- hitAnalysis measurement.txt --workPackageBasedAnalysis --save-cluster-binary --save-cluster-ascii --no-plots --no-save

4. Loading Data

The possible filters or variable parameters are defined in the "settings_measurement.ini" file, located in the same folder. The filters that are activated are for minimum and maximum photon count and a filter of energy in some of the plots. A copy of the used file is presented in the Section Settings Measurement.ini.

In order to work in Python with the generated flood maps, we need to extract the data, namely the four different calculated positions ("x" and "y" coordinates) and a flag indicating whether they are valid or not, from the output file ".DebugSingles". The condition of valid is set in the hitAnalysis software and currently corresponds to the fulfilment of the MainDie Neighbour condition. 000: MainDie is active (all clusters) and information only from MainDie, 100: MainDie and the Die next to it horizontally (depending on which pixel was activated), 010: MainDie and vertically, 111: MainDie, horiz, vertical and diagonal.

This data has be processed and ordered so it can be easily treated by the created routines in python.

Therefore, a new routine is developed for this purpose. This routine must be efficient and fast. The ".DebugSingles" files contain more than 100 million of clusters, so to reduce the use of RAM and be safer, each file is chunked in order to parallelize the process. The routine is saved in "ReadDebugSingles.py" and one only needs to indicate which path directory the .DebugSingles to analyze have and the path to the folder where the python-objects will be saved.

NOTE: The last cluster (event) of the .DebugSingles file is always corrupted and it does not have all the values, so this line is removed automatically.

After that, the corresponding hdf5 files will be created containing the information that we need for creating the final files.

4. Loading Data

A new framework has been created in python to process these data files. In the future, it can be fully adapted to MonoCal, but it runs independently now. The main routine is called 'NewHypmed.py'. The data is loaded using a specific function and then, the clusters are filtered by the validation flag and split into regions of interest (ROI). For every ROI a flood map (2d-histogram) is obtained and then saved with the library "pickle" which allows us to write on a file the python object.

The sequence of commands to obtain the final files to be used for the CIA is now given:

- python
- from NewHypmed import *
- Hyp = Hypmed()
- Hyp.sknet.load_raw_data_Pickles("PATHTODIRECTORYREAD")
- Hyp.valid_events()

- `Hyp.region_ana()`
- `Hyp.save_RefSection_into_file("SAVEDIRECTORY")`
- `Hyp.save_hist_into_file("SAVEDIRECTORY")`

5. Analysing Data

5.1. Peak Finder

Once we have the histograms saved on pickle files, we can now use a routine to find the peaks on the different histograms. It is based on the ROOT routine "SearchHighRes" and for every COG algorithm the parameters are changed (sigma, threshold, background and iterations of the convolution). The chosen values are shown in the following:

- `sigma = [3.3, 2.5, 1.5, 1.5] -> [000,010,100,111]`
- `threshold = [2, 6, 7, 4] -> [000,010,100,111]`
- `rmBackground = True -> forall`
- `convIter = 200 -> forall` with more tends to regroup peaks and with less some are missing
- `markov = False -> forall` smoothing function to remark peaks, but works for us, peaks with different heights and sigmas
- `mIter = 3 -> different number of iterations have been tried`

The output is an array of the "x" and "y" coordinates from the found peaks. Those coordinates are then ordered according to the "y" axis.

5.2. First Part of CIA

The first 20 peaks (clusters) are taken to be ordered by "x" axis. The closest to the left top corner is chosen as reference for the row. The next one must be within a range (a different one for every COG position algorithm): "y" position of previous peak \pm 0.1-0.2 mm. If not, the peak is saved for the next iteration. A scheme of this iteration is shown in Figure 1. This process is repeated 20 times because no more than 20 rows per ROI are expected. In case a peak or more could not be classified, they are saved for the second stage of the CIA in a dictionary.

For every histogram a dictionary of rows is created. An example of assigning rows in a ROI is shown in Figure 2.

5. Analysing Data




Figure 1: Scheme of iteration for obtaining rows.




Figure 2: Classification of rows from one ROI of 111-COG.

5.2.1. Assignment of Label to Crystals

In order to identify the peaks of the different ROIs with the corresponding crystal, a reference is needed. That means, it is necessary to assign the label of one crystal belonging to that ROI to, at least, one peak. With this first assignment, the rest of the peaks can be identified as we know the distribution of the crystals.

For that, the crystals must be identified with a label. As the orientation of the crystal will be consistent due to the disposition of the whole system (SiPM connections, Ring detector, etc.), the crystals will always have the same orientation on our measurements. The expected orientation can be observed in Figure 3.

Thanks to that fact, a consistent criteria for the labels is chosen from now on. The criteria follows a strictly order disposed by the position of the crystals. There are 3425 crystals and, for the sake of simplicity, the labels are numerical, so the labels will go from '0' up to '3424'. Using Figure 3 as a guide, the labels are set as follows.

The label '0' is given to the crystal at the top left corner corresponding to layer '1'. The next labels are assigned following the order of the row from left to right. After all the crystals of the first row belonging to the first and second layer are over, the next

5. Analysing Data

row belonging to the third layer follows. The assignment of the labels finishes with the last crystal at the bottom right corner being labelled as '3424'.

Apart from the label, a LUT_{labels} is created to order the labels into columns and rows. All the labels share the same row index, from '0' at the top until '70' at the bottom, which gives 71 rows in total.

However, the column index is different for those label belonging to the third layer respect to the other two layers. For both cases, the column index starts at '0' at the left side. For the third layer, the column index goes up to '30', as there are 31 columns in total for this layer. For layers 1 and 2, the column index goes up to '64', as there are 65 columns in total counting both layers.

A relation between the column indexes is given in Equation 1.

$$i_{12} = i_3 \cdot 2 + 2, \quad (1)$$

where i_{12} corresponds to the column index of the labels from the layers 1 and 2, and i_3 to the column index of the labels from the third layer.

5.2.2. First Identification

The chosen reference crystal per ROI for our algorithm is the one corresponding to the closest peak of the center of the ROI. This peak is assigned as it is the easiest to identify (less distortion and most stable area of ROI).

Since the crystal distribution is constant, a LUT_{ref} can be created with the labels corresponding to those crystals at the center of the ROIs for the different algorithms.

It is only necessary to look once at the ROI of the flood map for each algorithm to identify one of those crystal and the rest will be distributed systematically. By comparing the flood map with the expected distribution on Figure 3, one can find the peak which will be at the center of the ROI for each algorithm.

It is important to consider how the different algorithms work. In the case of the algorithm corresponding to '000' the ROI is concentrated on the center of the DPC. That means, that the peak at the center of the DPC will not suffer much distortion by the condition of the algorithm and it will correspond to the actual crystal that is over the center of the DPC. With the rest of the algorithms, it occurs similarly. We can observe this pattern on Figure ??.

Include image of 000,100,010,111 on one column showing a roi and the same are from expected floodmap!!!!!!!!!!!!!!

As a note: all the reference crystals belong to the third layer.

Thus, by identifying one crystal on one ROI per algorithm, we can add the rest of reference labels automatically to a LUT_{ref} .

5.2.3. Classification of Crystal

After this identification, the rest of the peaks belonging to the same row are classified. By knowing in which position of the row they are (as mentioned before, the row is sorted), by iterating over the column index and

5. Analysing Data

For example, if the peak in the center has the label 835, the label for the peak next to it (increasing the value of "x") in the same row would be 836. On the other direction, it would be 834. Then, the closest peak to this reference from the next row considering the "x" coordinate would be the reference label for the corresponding row. Here, we select how many rows and columns are accepted as valid pro ROI depending on the flood maps.

WRITE conditions and example

This information is saved in the general dictionary for the corresponding COG, where all the labels are variables of the dictionary. In case one peak is assigned to one label, this information is saved to the dictionary together with the layer to which belongs. In case there was already a peak from a different ROI, the new peak is appended. This case is known as ghosting, a peak which appears in two consecutive ROIs due to the fact that the crystal is between the regions.

To distinguish the peaks from the different layers, we can also use the reference as a guide. We know to which layer it belongs to and we know that if it is the third layer, the next one would be also belonging to the third layer and the reference for the next row would be a peak from the second layer. Next to the peak from the second layer, we have a peak from the first layer and so on. The general expected distribution of the peaks are shown in Figure 4.

Once we have assigned all the peaks from the dictionary of rows, the next ROI is iterated. After all the ROIs from the corresponding COG are obtained, a file is created with the dictionary including all the labels. A flow chart representing this procedure is shown in Figure 5.

5. Analysing Data




Figure 3: Expected distribution of centers of crystals for the different layers of the scintillator. The grey lines indicate the positions of the Pixels from the detectors. The dashed lines show the edges of the crystal. The pointed lines show the dimensions of the single crystals.

5. Analysing Data




Figure 4: Expected distribution of peaks from different layers.




Figure 5: Flow chart of the first part of the CIA.

5.3. Second Part of CIA

After the first classification, some of the peaks could be not correctly classified (as we can see in Figure 2) or directly not classified. That is why, a second check is done. Firstly, all the classified peaks from all the ROIs are taken and the median of the column from those valid ones belonging to the same column is calculated. In case they are within a specific range: median of column ± 0.2 - 0.3 mm (depending on which COG we are), they are confirm as valid. When that is not the case, the peak is marked as invalid.

The same procedure is followed with the medians of the rows. Now, all those peaks that are not valid, are then assigned to the closest column and row. If the new label is empty, the peak is assigned to it. If the label is invalid, we double check again the corresponding old and new peaks. Individually, each is checked whether they are within the specified range and then marked as valid.

At last, those peaks that did not have a label are also analysed. Their closest row and column are calculated and the same check is done. Finally, the new dictionary with the corrected labels is saved again in a file.

Now, in Figure 6 we can observe again the same ROI as shown in Figure 2, but with the correct classification of all the peaks.




Figure 6: Final classification of peaks with columns and rows. The orange dots are valid peaks and the blue dots are labels with two peaks.

5. Analysing Data

An example for a dictionary element follows:

```
dic_crystal[1715] =
{'center': 12: [[3.359644, 0.001895]], 13: [[4.560458, 0.00071]], (12 and 13 are ROI's id)
'id': 1715,
'layer': 3,
'row': 35,
'valid': True}
```

A simplified flow chart of this part is shown Figure 6.




Figure 7: Flow chart of the second part of the CIA.

6. Cluster Identification

Once the four COGs algorithms are passed to the CIA and four different dictionaries of labels are created, a cluster identification can be carried out.

To do so faster, a look up table (LUT) is created for all of them. A grid with a resolution of 0.1 mm is created for each COG algorithm. Each coordinate is then analysed to look for the closest label from the dictionary. The two closest peaks are saved in the LUT together with their relatives distances. If the distance is larger than 1 mm, the coordinate is not assigned to improve the classification later.

Example of lut element:

```
lud111[0,0] =
{'2CLOP': 'center': 12: [[-0.160166, -0.957595]],
'id': 1760,
'layer': 2,
'row': 36,
'valid': True,
'CLOP': 'center': 12: [[-0.160166, 0.001895]],
'id': 1712,
'layer': 3,
'row': 35,
'valid': True,
'QF': 0.026496736254681918}
```

Now, every single cluster can be round to 0.1 mm of resolution. These coordinates are looked in the four different LUTs from the COGs. If they have a label assigned and it is valid, the label is saved for comparison. When the valid labels agreed, the label is assigned to the cluster. When not, a coincidence between two or more of them has the priority to assign the corresponding label to the cluster. If none of them agreed, the quality factor (QF) extracted from the relative distances of the closest peaks has the priority. The quality factor is defined as:

$$QF = \frac{\text{lowest distance}}{\text{lowest distance} + \text{second lowest}} \quad (2)$$

The lower the QF, the better the peak is, because it means the relative distance to the first peak is smaller than the relative distance to the second peak and, that means, that in this COG this cluster is clearly identified to one peak. When the distances are similar, it means that the cluster is not clearly identified with only one peak. Thus, it is a worse identification.

A test of the cluster classification has been carried with 75 million clusters. The number of clusters assigned to each label is represented in Figure 8. Some of the crystals have more clusters assigned as others. On the edges it is normal to have more clusters assigned to the crystals because as we can observe some of the first and last rows are missing. They are not identified because they are collapsed in the flood maps and

6. Cluster Identification

therefore, they are not distinguishable. The other differences encountered could be due to the different filters applied. A deeper look into the data is needed to check whether there are problems in the classification.




Figure 8: Number of clusters in each label in logarithmic scale. Those with 0 clusters are not represented (first and last rows).

A. Settings Measurement.ini

```

1 #[frameCounter]
2 #bitWidth=8 #we use the 32 bit mec framecounter and thus do not need to
   correct the 8bit pdpc framecountner anymore

4 [dieConfiguration]
5 inhibitFraction=0.1

7 [clock]
8 frequency_MHz=195

10 [recPars]
11 randomizeSinglesInCrystalBin=false

13 #a lot of the following parameters can be optimized by using a special
   result extraction method that produces a debug print for the memory
   used by each workpackage (commandline options "--
   workPackageBasedAnalysis --analyseWorkPackageSize")
14 #this is the number of frames that are put into one workpackage (this
   should be 100 if raw data processed by the daps is used as the daps
   puts it also in 100 frame junks)
15 workpackageSize=50
16 #the number of reading threads is equal to the number of inout files (for
   MEC scanner it is usually 10)
17 #the following memory pool sizes are for each reading thread
18 #this is the size (byte) of one memory pool storing msg objects (without
   the message content)
19 readingThread_memory_msgs_poolSize=4000000
20 #this is the maximum number of memory pools for messages a single reading
   thread is allowed to allocate (if reading is faster than processing it
   will allocate this amount)
21 readingThread_memory_msgs_maxNumberPools=5
22 #the following two values are again for reading threads but this time for
   the memory used to store message content
23 readingThread_memory_msgMemory_poolSize=10000000
24 readingThread_memory_msgMemory_maxNumberPools=4
25 #this is the number of processing threads that are started
26 #processing threads will work on workpackages filled by the reading
   threads with the raw messages and will go through the complete data
   processing
27 #the number of threads can be chosen from 1 up to any number (it makes
   sense to keep it a bit smaller than the number of available cores [
   reading threads and result extraction thread(s)])
28 numberOfProcessingThreads=6
29 #the following values define the memory pool size and maximum number of
   memory pools per processing thread for hit, cluster and coincidence
   objects
30 processingThread_memory_hits_poolSize=8000000
31 processingThread_memory_hits_maxNumberPools=5
32 processingThread_memory_clusters_poolSize=32000000
33 processingThread_memory_clusters_maxNumberPools=5

```

A. Settings Measurement.ini

```

34 | processingThread_memory_coincidences_poolSize=1000000
35 | processingThread_memory_coincidences_maxNumberPools=5
36 | #result extraction (e.g. listmode files and plotting) is coordinated by a
   |   single thread that work on completely processed workpackages.
37 | #this thread will schedule plotting jobs for seperate threads to multi
   |   thread result extraction, this is not configurable atm

41 | #number of hits to be processed if the --max-hits option is set (is not
   |   100% accurate)
42 | NumberHitsToProcess=5e9
43 | #maximum number of ml poisson iterations
44 | MLPoissonMaxIterations=5
45 | #cluster window / ps
46 | clusterWindow=40000
47 | #coincidence window / ps
48 | coincidenceWindow=2000
49 | #method for the anger algorithm. possible values: allChannels,
   |   mainPixAndNeighborPixels, mainPixAndNeighborDies/onlyNeighbors,
   |   circleAroundMainHit, sameQuadrant, mainPixAndDirectNeighborPixels,
   |   mainPixAndAutoNeighborPixels
50 | angerMethod=mainPixAndNeighborPixels
51 | #if anger mehtod is circleAroundMainHit then this defines the radius
   |   around the main hit:
52 | angerMethod\angerRadius=10.
53 | #use the following method to identify the crystal: angerMap (default, if
   |   ommited), oneToOneCoupling
54 | crystalIdentificationMethod=angerMap
55 | #use a simple logarithmic saturation correction
56 | photonCountCorrection=saturation_log
57 | #apply the pixel gain if it has been calibrated
58 | usePixelGain=true
59 | #photon counts used for energy calculation. possible values: allChannels,
   |   mainAndNeighborPixels/onlyNeighbors, mainAndAutoNeighborDiesSum,
   |   mainAndNeighborDies, mainAndDirectNeighborDies, mainDie, mainPix,
   |   mlPhotonSum
60 | energyCalculationMethod=mainAndNeighborDies

62 | #cluster filter definitions: you can define cluster filters at different
   |   stages of the processing chain
63 | # -there are cluster filters used only for plots which do not influence
   |   the processing at all
64 | #   these filters are bound to a certain plot in the main.cc
65 | #   /clusterPlots - the standard filter used for clusterPlots using a
   |   filter
66 | #   /clusterPlots/[subIndex] - atm only index 0 is used for a filter
   |   without an anergy window
67 | # -/beforeCrystalIdentification - can be used for a fast photon cut (no
   |   neighbor information!!!! or energy is available !!!!!)
68 | # -/clusterFile - is only used directly before piping clusters into the
   |   cluster binary file per tile, the beforeCrystalIdentification has

```


A. Settings Measurement.ini

```
already been applied
69 # -/beforeCoincSearch - is applied before searching for coincidences,
    this may influence the random rate, if you filter before the coinc
    search as you may have a triple which loses one cluster
70 # -/coincSearch - is applied to coincident clusters, if no filter is
    applied before the coinc search one can reject all multi events
71 #
72 #possible filter settings (the /beforeCrystalIdentification does only
    support the photon cuts !!!!):
73 # name: string with name
74 # minPhotons, maxPhotons: integer values
75 # minEnergy, maxEnergy: floating point values
76 # neighborFilter: filter on the presence of the neighboring pixels of the
    main pixel; string value: none, directNeighbors, diagonalNeighbors,
    allNeighbors
77 # minMainPixRatio: the minimal required ratio of the photon value in the
    main pixel to the total number of photons in a cluster; floating point
    value [0., 1.]
78 # maxAngerSigma: maximal distance of a cluster to the nearest crystal
    anger position in terms of sigma; floating point value e.g. 3.

80 #clusterPlots: these are used for cluster plots and do not influence the
    processing, the mapping to plots is hardcoded in the main.cc!!
81 clusterFilter/clusterPlots/name="standard"
82 clusterFilter/clusterPlots/minPhotons=500
83 clusterFilter/clusterPlots/maxPhotons=5000
84 #clusterFilter/clusterPlots/minEnergy=411
85 #clusterFilter/clusterPlots/maxEnergy=561
86 clusterFilter/clusterPlots/neighborFilter=allNeighbors
87 clusterFilter/clusterPlots/0/name="noEnergy"
88 clusterFilter/clusterPlots/0/minPhotons=500
89 clusterFilter/clusterPlots/0/maxPhotons=5000
90 clusterFilter/clusterPlots/0/neighborFilter=allNeighbors
91 clusterFilter/clusterPlots/1/name="noEnergy_direct"
92 clusterFilter/clusterPlots/1/minPhotons=500
93 clusterFilter/clusterPlots/1/maxPhotons=5000
94 clusterFilter/clusterPlots/1/neighborFilter=directNeighbors
95 clusterFilter/clusterPlots/2/name="noEnergy_diagonal"
96 clusterFilter/clusterPlots/2/minPhotons=500
97 clusterFilter/clusterPlots/2/maxPhotons=5000
98 clusterFilter/clusterPlots/2/neighborFilter=diagonalNeighbors
99 #beforeCrystalIdentification
100 clusterFilter/beforeCrystalIdentification/name="
    beforeCrystalIdentification"
101 #clusterFilter/beforeCrystalIdentification/minPhotons=1000
102 #clusterFile
103 clusterFilter/clusterFile/name="clusterFile"
104 clusterFilter/clusterFile/neighborFilter=allNeighbors
105 clusterFilter/clusterFile/minPhotons=500
106 clusterFilter/clusterFile/maxPhotons=5000
107 #clusterFilter/clusterFile/minEnergy=
108 #clusterFilter/clusterFile/maxEnergy=
```

A. Settings Measurement.ini

```
109 #beforeCoincSearch: this filter is applied BEFORE the clusters are
    processed by the coincidence finder
110 clusterFilter/beforeCoincSearch/name="beforeCoincSearch"
111 clusterFilter/beforeCoincSearch/minPhotons=500
112 clusterFilter/beforeCoincSearch/maxPhotons=5000
113 #clusterFilter/beforeCoincSearch/minEnergy=411
114 #clusterFilter/beforeCoincSearch/maxEnergy=561
115 clusterFilter/beforeCoincSearch/neighborFilter=allNeighbors
116 #filter for coincidences: minTiles, maxTiles, minLORLength (/ mm),
    minTileDistInPhi AND/OR minSPUDistInPhi (USE ONLY FOR MEC SCANNER WITH
    10 SPUS!!!!) (specifies the minimum number of tiles/SPUs in BETWEEN
    the two coincident tiles/SPUs)
117 coincidenceFilter/minTiles=2
118 coincidenceFilter/maxTiles=2
119 #coincidenceFilter/minLORLength=500
120 #coincSearch: this filter is applied after the coincidence search to every
    cluster of a valid coincidence.
121 # You only need to implement stricter filters than the filter
    before search, because all clusters must have already passed the
    previously applied cluster filter
122 clusterFilter/coincSearch/name="afterCoincSearch"

124 [sinogramPars]
125 # Half the crystal width thanks to intercept theorem:
126 # The distance between the two LORs L_0,300 and L_0,301 is 0.5 mm
127 # at half of the length of the LORs for a crystal width of 1 mm
128 displacementResolution=2
129 # Inner diameter
130 displacementSize=217.6
131 # sin(pitch/radius) * 180/pi
132 # radius instead diameter to prevent checker board pattern
133 angularResolution=2.1
134 NRings=8
135 # These are the upper edges of the axial crystal bins
136 ringsLUT=-12, -8, -4, 0, 4, 8, 12, 16
```